

Required materials

- Coating buffer: 0.1 M Sodium Carbonate, pH 9.5 (8.4 g NaHCO₃, 3.56 g Na₂CO₃, 1 litre de-ionized water).
- Wash-buffer: PBS, 0.05% Tween20.
- Blocking-buffer/Assay diluent: PBS, 1% BSA
- Capture antibody: Purified anti-mouse IgE, IgG1, or IgG2a antibodies.
 - Anti-mouse IgE (Citeq part number, HE.02.01), dissolve in 0.5 ml coating buffer, dilute 20x.
 - Anti-mouse IgG1 (Citeq part number, HE.02.02), dissolve in 0.5 ml coating buffer, dilute 20x.
 - Anti-mouse IgG2a (Citeq part number, HE.02.03), dissolve in 0.5 ml coating buffer, dilute 20x.
- Optional: Reference mouse serum (for relative quantification or as a positive control)
- Recommended: Non-immunized mouse serum (for use as a negative control)
- Biotinylated Der p 1, Der p 2, D. pteronyssinus extract or D. farinae extract (available from Citeq):
 - D. pteronyssinus extract Biotinylated (Citeq part number, 02.01.88), dissolve in 0.5 ml Buffer (recommended PBS), dilute 50x.
 - D. farinae extract biotinylated (Citeq part number, 02.02.88), dissolve in 0.5 ml Buffer (recommended PBS), dilute 50x.
 - Der p 1 biotinylated, (Citeq part number, 02.01.73), dissolve in 0.5 ml Buffer (recommended PBS), dilute 20x.
 - Der p 2 biotinylated, (Citeq part number, 02.01.74), dissolve in 0.5 ml Buffer (recommended PBS), dilute 50x.
- Streptavidin-horseradish peroxidase conjugate.
- TMB substrate solution
- Stop solution: 2.5 M H₂SO₄.
- Flatbottom 96 Microwell plates.

Required equipment

- Microplate reader with software
- Microplate washer (optional)

Method

1. Coat micro wells with 100 µl per well anti-IgE (or IgG1/ IgG2a) capture antibody diluted in coating buffer. Incubate overnight at 4°C.
2. Aspirate wells and wash 3 times with 300 µl/well wash-buffer.
3. Block plates with 300 µl blocking-buffer. Incubate at room temperature for (at least) 1 hour.
4. Aspirate wells and wash 1 time with 300 µl/well wash-buffer.
5. Dilute serum samples and controls/references in assay diluent. Pipette 100 µl of each in duplicate into the wells. Incubate for 2.5 hours at room temperature.
6. Aspirate wells and wash 5 times with 300 µl/well wash-buffer.
7. Add 100 µl/well of diluted biotinylated HDM/Der p 1/Der p 2 to each sample and incubate for 1 hour at room temperature.
8. Aspirate wells and wash 5 times with 300 µl/well wash-buffer.
9. Add 100 µl of diluted Streptavidin-Peroxidase to each well, incubate for 30 minutes at room temperature.
10. Aspirate wells and wash 7 times with 300 µl/well wash-buffer.
11. Add 100 µl of substrate solution to each well.
Incubate plate at room temperature in the dark until:
 - a. blanks are staining or
 - b. until the positive control stains clearly.
12. Add 50 µl of stop solution to each well.
13. Read absorbance at 450 nm within 30 minutes of stopping reaction. Use A 650 as a reference filter.